

CROWN

Operation Manual



CROWN

Thank you for your purchase our CROWN head massager.
Before operation the CROWN head massager, please read
the User's Manual thoroughly for the correct methods of usage,
especially the notes on safely. This manual should be kept
available for future reference after you have read it through

CONTENTS

Names of Functions of components	1
Preparations	2
Function & Specialty	3
NOT Suitable Users	4
Storage & Maintenance, Warning	5
Important Safety Instructions	6
Specifications	6

Suitable Users

Suitable Users

- Those who are under medical care or are feeling unwell
- Those who are suffering from malignant tumors
- Those who are suffering from heart diseases
- Women who are pregnant
- Those who are under medical rest or rest-cure ordered
By doctor
- Those who have high fever

Warning


- Close supervision is necessary when this appliance is used by, on or near children, invalids, or disabled persons
- Do not using who are head & hair is wet .
- Do not still to using who feeling is unwell
- Persons who suffer from brain diseases or have done brain operations should not use the massager.
- Take off the contact lens ,hairpin and hair decoration before using the massager.
- Do not clean the device with water, If it needs cleaning, wipe it with aseptic cotton dipping alcohol.
- Do not hold the cable of the device only when you use the massager. Do not disassemble the device yourself.
- Please keep away from direct sunlight or heating appliances
Like heat stove
- Do not walk outside or drive when you use the massager.
- Switch off the power when completed the massage

⚠ Important Safety Instructions

- Inspect the power supply. Do not turn on the device if the power voltage exceeds the specified range.
- Inspect eye patch and cable regularly, If any damage or breaking occurs, stop using the massager immediately
And send it to our after service center.

Names of Function of Components


Function & Specialty



- Ⓐ Sliding Size Extend
- Ⓑ Air Cord
- Ⓒ Head Vibration.
- Ⓓ Massage Main Unit

Names of Function of Components

Function & Specialty



Program Reset

The reset mode is minutes with air pressure. The massage mode work repeatedly.

MUSIC

Press the MUSIC selection button to change the music.

AIR PRESSURE switch

Automatic running mode of air pressure massage.

TIME switch

Press the TIME selection button to change the working time (5.10.15 minutes)

VIBRATION switch

Automatic running mode of vibration.

INFRARED

Press the INFRARED button to select the infrared On /Off.

Storage & Maintenance

Storage

- Make sure the power is turned off and pull out the plug from the outlet when not using the device..
- Use dry and soft cloth to clean up the cushion
- Do not use diluents, volatile oil or bleaching powder will cause Casing Change color or damaged
- Do not use Lubricant on this product
- Keep the device where it is dry and there's less dust

Product Specification

Model Name	: AIR CUSHION HEAD MASSAGER
Model no	
Frequency	: DC 6V
Power consumption	: 6W
Rating Time	: 15 minutes

CROWN